


# Filter Funnel



**Product Categories:** [Funnels Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/filter-funnel/>

## Product Description

Filter Funnel - Laboratory Glassware Equipment

A filter funnel is a laboratory funnel used for separating solids from liquids via the laboratory process of filtering. In order to achieve this, a cone-like shaped piece of filter paper is usually folded into a cone and placed within the funnel.

The suspension of solid and liquid is then poured through the funnel.

**Features:**

High strength

Precise Design

Easy to handle

Dia (cms) 5, 6.25, 7.5, 10.15

Accurate measurement

Transparent