

Filter Paper



Product Categories: [General Chemistry Labware](#)

Product Tags: [chemistry lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/filter-paper/>

Product Description

Filter Paper - Chemistry Laboratory Equipment

The main aim of filter paper is to separate soil particles from a liquid purpose of filtration to clarify liquid. to recover solid. Filter paper is the only media that yields exceptional clarity all grades of filter paper are made of highest quality of raw materials with selective grade & that is why it is used for strength, filtration speed retentivity purity.

Specifications:

Low ash content

Non reactive with any organic compound

Varied sizes, microns and other related specs.

Ream 46x 57cm

Regular quality.

Superior quality

Filter paper is a semi-permeable paper barrier placed perpendicular to a liquid or air flow.

It is used to separate fine substances from liquids or air.

Available in size - 11cm, 12.5cm, 15cm