

## Filtration Assembly Borosilicate Glass



**Product Categories:** [Filtration Assembly Borosilicate Glass](#), [Laboratory Glassware](#)

**Product Tags:** [glassware](#), [Laboratory Equipment](#)

**Product Page:**

<https://www.labappara.com/product/filtration-assembly-borosilicate-glass/>

### Product Description

Filtration Assembly Borosilicate Glass - Laboratory Glassware Equipment

Complete, Filter Flask with Socket and buchner

Filtration Assembly Set Manufactured from Borosilicate glass 3.3 which is heat resistant. Particular suitable analysis of microscopic residue Particular control Borosilicate 3.3 sintered disk for filtration to be used with 47/50 mm dia. Membrane can be autoclave at 120o for 20 minutes.

### Features:

Easy to use

Long life

Excellent filtration

Crack proof and transparent

funnel with Sintered Disc of Porosity G-1 or G-2 or G-3 or G-4, having cone.

Cat. No.

Flask Cap.

Buchner Funnel

Joint

327/1

250

80

24/29

327/2

500

200

24/29

327/3

1000

500

29/32

327/4

2000

1000

29/32