


Frequency Modulation & Demodulation Trainer (Advanced Model)



Product Categories: [Advance Communication Lab](#), [Engineering Equipment](#)

Product Page: <https://www.labappara.com/product/frequency-modulation-demodulation-trainer-advanced-model/>

Product Description

Frequency Modulation & Demodulation Trainer (Advanced Model)

Features:

- Reliable
- Durable finish standards
- Excellent performance

Specifications:

- +/-12V/350mA DC Power Supply
- Carrier Generator 75 KHz
- Modulating signal Generator 10 KHz
- Frequency Modulator
- Frequency Demodulator
- Set of 2 mm Patch Chords & Manual