

Funnel-Powder, Short Plain Stem



Product Categories: [Funnels Borosilicate Glass](#)

Product Tags: [borosilicate glass](#), [laboratory glass ware](#), [laboratory glassware](#)

Product Page:

<https://www.labappara.com/product/funnel-powder-short-plain-stem/>

Product Description

Funnel Powder, Short Plain Stem- Laboratory Glassware Equipments

The power funnels are widely used for separating the debris from powder.

Made of Borosilicate 3.3 heat resistant glass

Powder Funnel with cone.

Features:

Easy to use

Durable

Perfect shape

Sturdy design

Accurate marking

Crack resistance

Smooth edge

Available In diameter- 55 mm, 80 mm, 100 mm, 125 mm, 150 mm