

# Galvanometer



**GALVANOMETER**

**Product Categories:** [Lab Meters](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/galvanometer/>

## Product Description

### **Galvanometer- Physics Equipments**

A galvanometer is an electromechanical instrument for detecting and measuring electric current.

The most common use of galvanometers was as analog measuring instruments, called ammeters, used to measure the direct current (flow of electric charge) through an electric circuit.

It is commonly used In Schools, Colleges and Universities Laboratories.

Dc mounted in beautiful Bakelite / acrylic stand with accuracy  $\pm 0\%$  F.S.D.

With lock Type terminal packed in card box.

All common ranges available.

Available Models – MO65, MR80, MR100