


## Glass Slab



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/glass-slab/>

### Product Description

#### Glass Slab - Physics Optical Equipment

Rectangular, all sides polished, made from slightly greenish glass free from bubbles. These are used for pin method experiments. Glass slab from Sheet Glass Bubble less, prepare from Belgium plate, highly polished sides, English glass.

#### **Specifications:**

Material : Glass

Glass Color : Transparent

Use : Lab use

#### **Features:**

Optimum finish

Crack resistance

High strength

**SIZE:-**

- (A) 75 x 50 x 12 mm
- (B) 75 x 50 x 18 mm
- (C) 100 x 60 x 25 mm