

Gynecological Delivery Bed



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/gynecological-delivery-bed/>

Product Description

Gynecological Delivery Bed - Hospital Equipment

Constructed from MS tubular steel material. Three section top with adjustable head section, seat section and leg section.

Top is upholstered and covered with water proof synthetic material.

Complete with SS bowl 32 cm dia.

Standard padded leg holder finish in epoxy power coated.

Overall Size - 1830 L X 550 W X 800 Hmm