

# H' Parameter Of PNP Transistor In Common Emitter Mode Apparatus



**Product Categories:** [Characteristics And Application](#), [Electronics](#)

**Product Page:** <https://www.labappara.com/product/h-parameter-of-pnp-transistor-in-common-emitter-mode-apparatus/>

## Product Description

H' Parameter Of PNP Transistor In Common Emitter Mode Apparatus

**Features :**

Instrument comprises of Two DC Regulated Power Supplies 0-10VDC/150mA, 0-1VDC/150mA, Four round meters, One PNP Transistor mounted behind the panel, connections of Supplies, Meters & Transistor brought out at 4mm Sockets.