

Hand and Table Microtome



Product Categories: [Biology Lab](#), [Educational equipments](#), [Microtome](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/hand-table-microtome/>

Product Description

Hand and Table Microtome - Biology Laboratory Instruments

Hand and Table Microtome for Biology Lab - Hand and Table Microtome

- Microtomes Biology Stage dia. 80mm. vertical lift 20 mm. and least count 0.01 mm,
- Can be securely attached to the edge of bench, table or work surface with the help of fixing clamp.

Features -

- With table clamping device and Plano razor.
- Feed adjustment in step of 10 microns.
- Optimum performance
- Less maintenance