

# Hospital Examination Couch Technical Specification



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:**

<https://www.labappara.com/product/hospital-examination-couch-technical-specification/>

## Product Description

Hospital Examination Couch Technical Specification - Hospital Equipment

Constructed from CRCA steel sheet.

Manually Adjustable Head Section.

Two individual cabinets with attached drawers

Epoxy power Coated.

**Overall Size - 1830 L X 60 W X 810 H**