

## Hospital Full Flower Bed



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/hospital-full-flower-bed/>

### Product Description

Hospital Full Flower Bed - Hospital Furniture & Equipment

Bed's frame work of strong and precise steel tubes having 'H' type legs on 10 cm castor two with brakes.

Having four section top of perforated CRCA sheet duly machine press bent with mattress guards.

Backrest and knee rest section maneuvered by two separate screws from foot end.

ABS panels Board.

IV Rod locations.

**Overall Size** – 208L X 90W X 60H cm.