

## Hospital Warden Bed with Back Rest



**Product Categories:** [Hospital Equipment & Furniture](#)

**Product Tags:** [hospital](#), [medical](#)

**Product Page:** <https://www.labappara.com/product/hospital-warden-bed-back-rest/>

### Product Description

Hospital Warden Bed with Back Rest - Hospital Equipment & Furniture

Strong rectangular steel tube frame work, perforated CRCA sheet top bed.

Adjustable brackrest upto 80 deg.

Fix lever handle.

Tubular head and leg bows of unequal height with vertical tube support.

Power coated finishing.

**Overall Size** – 208L X 92W X 60H