


Induced Current Apparatus



Product Categories: [Educational equipments](#), [Electromagnetism](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/induced-current-apparatus/>

Product Description

Induced Current Apparatus - Physics Equipment

Induced Current Apparatus has the bulb in secondary and break arrangement in primary to show the phenomenon of self and mutual induction and principle of induction coil, workable on 6 volts.

A soft iron core fits the cylinder of the primary.

Complete with terminals.