

## Ionization Potential Of Mercury Using Thyatron Valve



**Product Categories:** [E/M Apparatus & Valve Experiments](#), [E/M Apparatus & Valve Experiments](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:**

<https://www.labappara.com/product/ionization-potential-of-mercury-using-thyatron-valve/>

### Product Description

#### Ionization Potential Of Mercury Using Thyatron Valve

**Features:-**

Instrument comprises of DC regulated power supply 0-15VDC/150mA & AC power supply 6.3 VAC/1Amp, two round meters for voltage & current measurement, valve diagram is printed on front panel, connection of supplies and valve brought out at 4mm sockets.