

Junior Medical Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/junior-medical-microscope/>

Product Description

Junior Medical Microscope - Scientific Equipment Magnification 20 x to 1500 x

Stand Robust and stable, inclinable through 90 degrees. Mechanical tube length 160 mm.

Focusing Separate knobs for coarse and fine motion.

Stage Fixed square stage 120 x 120 mm size with mechanical stage having travel range 30 x 76 mm with graduated scales and vernier.

Illumination Plano concave mirror on gable mount. Optionally available attachable 20 w lamp for 220 v or 6 v, 20 w with transformer.

Nose piece Quadruple- revolving

Objectives 10 x & 45 x (s. L.) with 100 x oil immersion.

Eye pieces Huygenian 5 x or 15 x, wf 10 x.