

## L.D.R Characteristic Apparatus



**Product Categories:** [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:** <https://www.labappara.com/product/l-d-r-characteristic-apparatus/>

### Product Description

L.D.R Characteristic Apparatus - Electronics Laboratory Equipment

L.D.R Characteristic Apparatus with build in regulated power supply.

With build in regulated power supply and with 2 meters.

**Objective :** To study Characteristics of LDR

#### **Specification:**

Power Supply : 0-10 Volt DC

Voltmeter : 0-10 Volt DC

Ammeter : 50 mA

Lamp & Wood Plank : 100W , 50cm Wooden Plank

Size : 11"x9"x3" inch

#### **Features :**

Excellent durability

Smooth functioning

Low maintenance cost