

Labour Bed



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/labour-bed/>

Product Description

Labour Bed - Hospital Equipment & Furniture

Constructed with MS rectangular tube steel material. Adjustable head and leg section on castors, with two breaks. Top is upholstered and covered with water proof synthetic material.

Complete with SS bowl 32 cm dia.

Standard padded leg holder

Finish in epoxy power coated

Overall Size - 1830 L X 610 W X 800 H mm