

# LED Characteristics Apparatus



**Product Categories:** [Characteristics And Application](#), [Characteristics And Application](#), [Electronics](#), [Engineering Equipment](#)

**Product Page:** <https://www.labappara.com/product/led-characteristics-apparatus/>

## Product Description

LED Characteristics Apparatus

**Features :**

Instrument comprises of DC Regulated Power Supply 0-3 VDC/150mA, two round meters for voltage & current measurement, LED mounted on the panel, connections of Supplies & LED brought out at 4mm Sockets.