

## Magnifiers ( Imported)



**Product Categories:** [Educational equipments](#), [Light & Optics](#), [Magnifier](#), [Optical Instruments](#), [Physics Equipments](#)

**Product Tags:** [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/magnifiers-imported/>

### Product Description

Magnifiers ( Imported) - Physics Optical Equipment

The Reading Magnifier also can be used for outdoor activities. 2-in-1 lens mineral glass Magnifying at 3X, 6X magnification power to help you read the smallest lettering on maps and menus.

The hand-held magnifier design involves a handle that is attached to a lens through which the user will view an object - whether it is text or an image Hand-held magnifiers are the one of the simplest type of magnifiers to use and are also the most familiar. Nearly everyone has used a hand-held magnifier at some point in time.

<b>Diameter</b>	<b>Focal length</b>
-----------------	---------------------

75mm	20-100cm
------	----------

100mm	20-100cm
-------	----------