


MCT Box



Product Categories: [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/mct-box/>

Product Description

MCT Box - Scientific Laboratory Equipment

MCT Boxes are moulded in polycarbonate and are therefore strong as well as autoclavable. One size can hold MCT's of the other can hold MCT's . These boxes are designed to work at temperatures ranging from . The transparent cover provides a good view of the contents inside & is numbered for better inventory control. These boxes are designed to work at temperatures ranging from -190 C to 12 C.

Material : Polycarbonate