

Medical Microscope



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [biology lab equipment](#)

Product Page: <https://www.labappara.com/product/medical-microscope/>

Product Description

Medical Microscope - Laboratory Equipment

Stand : Stable & Robust, traditional Horse Shoe base body. Body inclinable through 90°, Standard tube length of 160 mm

Focusing : Separate coarse and fine focusing knobs with graduation on one side. Coarse movement by Rack & Pinion and Fine movement by Screw Lever Mechanism

Nosepiece : Quadruple Revolving Nosepiece is provided with positive click stop for proper centering

Stage : Fixed Stage size 125 X 125 mm with detachable mechanical stage for X and movement of slide up to 25 X 75 mm with graduated scales and verniers

Illumination : illumination through a 50mm diameter planoconcave mirror mounted on gimbals

Objectives : Three Achromatic Objectives : 10 X 40 X (SL) & 100 X (SL) Oil

immersion

Eyepieces : Two Huygenian Eyepieces : 5 X & 10 X

Condensor : Sub Stage Abbe condenser with N.A. 1.25 moves up and down by rack & pinion. Irs diapharagm and filter holder are mounted with the condenser

Safety Device : Adjustable down movement stopper for preventing accidental damage