

Magenta Box



Product Categories: [Laboratory Plastic Ware](#), [Sorage Boxes Carriers](#)
Product Tags: [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)
Product Page: <https://www.labappara.com/product/megenta-box/>

Product Description

Magenta Box - Laboratory Equipment

Magenta Box:- Magenta Boxes very often used in tissue culture & agricultural research experiments. Main body of this Magenta Box is made of tough & transparent poly-carbonate which provides excellent view of the contents inside while the cover is made of Polypropylene. This autoclavable Magenta Box is square in shape giving itself a space saving features.

Material : Poly-carbonate / Polypropylene

Particulars	Packing
Magenta Box	12 Pcs