

## Melting Point Apparatus



**Product Categories:** [Laboratory Equipments](#), [Melting Point Apparatus](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#), [Manufacturer of Laboratory equipment](#), [math lab instruments](#), [physics lab instruments](#), [science lab equipment](#)

**Product Page:** <https://www.labappara.com/product/melting-point-apparatus/>

### Product Description

#### Melting Point Apparatus - Scientific Laboratory Equipment

A metallic heating block is provided with holes to receive three capillary tubes and a stem type thermometer. The body is finished in stoving enamel paint. The apparatus can measure melting point up to 3500 C. The block is heated by superior heating elements and temp. is controlled by a heating regulator. A shadow less lamp has been provided to illuminate the capillaries and focus-able magnifier facilitates in observations. The whole unit is housed in a aluminum cast casing and provided with a control panel. Provided with indicator lamps, heat energy regulator, on-off switch, cord-plug etc. the equipment is to work on 220 V AC 50 Hz single phase.

#### Features:

Easy installation

Determine melting points of single solids

Smooth finish