

## Monocot Stem Anatomy, TS And LS



**Product Categories:** [Educational equipments](#), [Educational Models](#), [Zoology Models & Botany Model](#)

**Product Tags:** [biology lab](#), [botany model](#), [fibre model](#), [science model](#), [zoology](#)

**Product Page:** <https://www.labappara.com/product/monocot-stem-anatomy-ts-ls/>

### Product Description

#### **Monocot Stem Anatomy, TS/LS- Scientific Educational Model**

Exhibits the various tissues and the scattered closed and collateral vascular bundles in transverse and longitudinal sections in maize. The large pitted vessels, spiral and annular vessels show the anatomy of monocot stem. Mounted on board with numbered key card.