

## Motorized Gyroscope Apparatus



**Product Categories:** [Engineering Equipment](#), [Theory Of Machine Lab](#)

**Product Page:**

<https://www.labappara.com/product/motorized-gyroscope-apparatus/>

### Product Description

#### Motorized Gyroscope Apparatus

This is to demonstrate the relationship between applied torque and the rate and direction of rotation. Motorized Gyroscope This is motorized and mounted on a frame such that free to rotate about three perpendicular axis. This facilitates to observe the laws of stability and justification of equation.

#### Technical Details:

Stop Watch

Electronic

Variable speed of standard make

Supplied with Speed Control Unit

Weights: 2 kg, 1 kg,  $\frac{1}{2}$  kg

Stainless Steel

Diameter: 300 mm x 12 mm thick precisely balanced which can be rotated in 3 mutually perpendicular axis

Motor

Accurately marked scale & pointer to measure precession rate

Material

Disc