Multiple Adapters, Two Necks at 45°



Product Categories: Adapters Borosilicate Glass, Laboratory Glassware

Product Tags: borosilicate glass, laboratory glassware

Product Page:

https://www.labappara.com/product/multiple-adapters-two-necks-45/

Product Description

Multiple Adapter With Two Neck at 45°- Laboratory Glassware

Multiple in with Three Necks, Two Parallel and one at 45°.

Available In Socket Size: 19/26

Available In Cone Size: 19/26, 24/29, 34/35