

# Narrow Mouth Bottle



**Product Categories:** [Bottles And Asssesories](#), [Laboratory Plastic Ware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/narrow-mouth-bottle/>

## Product Description

Narrow Mouth Bottle - Scientific Laboratory Equipment

This autoclavable bottle, made of polypropylene, can be used for storing & transferring chemicals and other reagents safe & sound.

Range in size from narrow mouthed to wide mouthed and often a glass or plastic funnel is needed to properly fill a reagent bottle from a larger or equal.