

Newton Color Disc



Product Categories: [Microscope](#), [Microscopes](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/newton-color-disc/>

Product Description

Newton Color Disc: Physics Lab Equipment

Newton Color Disc for demonstrating that white light may be composed from all the spectral colors by rotation of multi color disc.

Comprising of a disc 200 mm, mounted on metal stand and driven by belt with the help of driving wheel fitted with a handle, can be clamped on table.

Specifications:

Disc 200 mm

Complete assembly mounted on a stable base.