

NEWTON'S LAW OF COOLING



Product Categories: [Heat And Thermodynamics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/newtons-law-cooling/>

Product Description

Newton's Law of Cooling - Physics Equipment

Newton's Law of Cooling explains heat transfer between two different bodies. This concept forms the basis of thermodynamics, hence understanding this concept is critical. This apparatus goes a long way in accomplishing that as it is simple to use and accurate in its performance.

Consist of two copper calorimeters resting inside a trough which stand in an outer tank

The latter is fitted with a wooden top with two holes for thermometer.

Supplied without thermometer, but with stirrer.

Available in -

GI sheet

Brass Sheet