

Oak Ridge Centrifuge Tube



Product Categories: [Centrifuge Bottles And Tubes](#), [Laboratory Plastic Ware](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/oak-ridge-centrifuge-tube/>

Product Description

Oak Ridge Centrifuge Tube - Scientific Laboratory Plastic ware

Oak Ridge Centrifuge Tubes, made of Polypropylene, have excellent contact clarity and can be used upto 50000xG. The screw cap is fitted with silicon O-Ring making it totally leakproof. These tubes should not be used below 4°C.

Particulars	Packing
30 ml.	12 Pcs
50 ml.	12 Pcs
100 ml.	12 Pcs