


## PCR Tube Rack



**Product Categories:** [Clamps, Holders Racks Stands](#), [Laboratory Plastic Ware](#)

**Product Tags:** [biology lab equipment](#), [chemistry lab instruments](#)

**Product Page:** <https://www.labappara.com/product/pcr-tube-rack/>

### Product Description

PCR Tube Rack - Scientific Laboratory Equipment

This Rack, as the name suggest, is used for keeping PCR Tubes. Top & base of this rack are made of tough & transparent polycarbonate while the middle portion is made of polypropylene. The tube holding spaces are uniformly placed so as to facilitate the use of multi pipettes. A no. of racks can be sacked together to save shelf space.

Material : Polycarbonate / Polypropylene

Rack for 96 PCR Tubes of 0.2 ml.