

## Pipettes, Measuring



**Product Categories:** [Laboratory Glassware](#), [Pipettes](#) , [Volumetric Flasks And Measuring Cylinder](#)

**Product Tags:** [glassware](#), [Laboratory Equipment](#)

**Product Page:** <https://www.labappara.com/product/pipettes-measuring/>

### Product Description

Pipettes, Measuring - Laboratory Glassware Equipment

**"Labappara"**. Accuracy as per with works certificate. Color coded. (These pipettes are calibrated to deliver from any graduation line down to zero at the jet, with the last drop expelled by blowing). Pipette Measuring made of Borosilicate Glass.

#### Features:

Accurate measurement

Durability

Fine finishing

Easy to handle

Unique design

Cat. No.	Capacity	Sub. Div.	Tolerance
MI.	ml.	± ml.	
39/1	0.1	0.01	0.006

39/2	0.2	0.01	0.006
39/3	0.5	0.02	0.006
39/4	1	0.01	0.007
39/5	1	0.1	0.007
39/6	2	0.02	0.01
39/7	2	0.1	0.01
39/8	5	0.05	0.03
39/9	5	0.1	0.03
39/10	10	0.1	0.05