

## PN Junction Diode Characteristics Apparatus



**Product Categories:** [Characteristics And Application](#), [Electronics](#)

**Product Page:**

<https://www.labappara.com/product/pn-junction-diode-characteristics-apparatus/>

### Product Description

#### PN Junction Diode Characteristics Apparatus

##### Features:

Instrument comprises of the following:

- Two DC Regulated Power Supplies 0-3VDC/150mA & 0-30 VDC/150mA
- Two dual range round meters for voltage & current measurement with connections brought out at 4mm sockets
- PN Junction Diodes mounted behind the panel
- Connections of Supplies, Meters & Diodes brought out at 4 mm Sockets.
- Front Panel made of thick Bakelite sheet
- One ON/OFF Switch with jewel light is provided on the front panel.
- Made of Heavy duty metal box construction

Accessories Supplied:

- Required number of Patch Chords
- Instruction Manual