

Powder Funnel



Product Categories: [Funnels Borosilicate Glass](#)

Product Tags: [borosilicate glass](#), [laboratory glassware](#)

Product Page: <https://www.labappara.com/product/powder-funnel/>

Product Description

Powder Funnel- Laboratory Glassware Equipments

The power funnels are widely used for separating the debris from powder.

Made of Borosilicate 3.3 heat resistant glass

Powder Funnel with cone.

Features:

Easy to use

Durable

Perfect shape

Sturdy design

Accurate marking

Crack resistance

Smooth edge

Powder Funnel With 24/29 cone at stem-

Available in Cone Size- 100, 125, 150