


Precision Epidiascope



Product Categories: [Audio- Visuals](#), [Educational equipments](#)

Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)

Product Page: <https://www.labappara.com/product/precision-epediascope/>

Product Description

Precision Epidiascope - Scientific Laboratory Educational Equipment

An efficient convertible projector equipped with high grade anastigmatic triplet white lenses & surface coated reflecting mirror.

Feature -

Providing brilliant projection of opaque objects as well as colored slides.

Equipped with a noiseless cooling system.

Providing undisturbed projection for long duration.

For use as an epidiascope, a slide projector, a film strip projector.

Micro slide projector with the aid of different accessories along.

Complete with 225v-100w imported tube.