

Projection Slides



Product Categories: [Audio- Visuals](#), [Educational equipments](#)
Product Tags: [biology lab equipment](#), [chemistry lab instruments](#)
Product Page: <https://www.labappara.com/product/projection-slides/>

Product Description

Projection Slides - Scientific Laboratory Educational Equipment

A slide projector is a device that is used to view photographic slides by using optical and mechanical methods. It contains an electric light bulb.

Feature-

Mounted in plastic slides

With heavy transformer

Supplied in moulded luggage carrying case.

Focusing lenses.

Reflector and condensing lenses

Light passes through the transparent slide and lens, and the resulting image is enlarged and projected onto a screen.