

Reagent Bottle (wide mouth)



Product Categories: [Bottles And Asssesories](#), [Laboratory Plastic Ware](#)

Product Tags: [Laboratory Equipment](#), [Laboratory plastic ware](#)

Product Page: <https://www.labappara.com/product/reagent-bottle-wide-mouth/>

Product Description

Reagent Bottle (wide mouth) - Scientific Laboratory Equipment

Wide Mouth Reagent Bottle:- Wide Mouth Reagent Bottle made of Polypropylene. These have unique features which include Rigidity, Translucent and have Excellent Contact Clarity

Material : Low Density Polyethylene & Polypropylene

These Wide Mouth Reagent Bottles are available in Polypropylene as well as Polyethylene.

The PP bottles are rigid, translucent and have excellent contact clarity.

These bottles can be autoclaved filled as well as empty.

The Polyethylene bottles too have excellent strength and are unbreakable.

The Polypropylene cap has a built in seal ring to make these bottles leakproof.