Reduction Adapters



Product Categories: Adapters Borosilicate Glass, Laboratory Glassware

Product Tags: glassware

Product Page: https://www.labappara.com/product/reduction-adapters/

Product Description

Reduction Adapters - Laboratory Glassware Equipment

All reduction adapters are made from borosilicate glass and joints are interchangeable with other apparatus. An outer interchangeable joint at the top (Socket) and a smaller inner interchangeable joint (Cone) at the base.

Available in socket size - 14/23, 19/26, 24/29, 29/32, 34/35, 40/38

Features:

Chemical resistance Leakage proof Flawless finish