

Riders (5mg to 10mg)



Product Categories: [Measurements](#), [Mechanics](#), [Physics Equipments](#)

Product Tags: [chemistry lab instruments](#), [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/riders-5mg-10mg/>

Product Description

Riders (5mg to 10mg) - Physics, Chemistry Equipment

A Rider is a small piece of aluminum, from which various small weights can be attached to a balance to determine the weight of small, light objects.

Features:

Accuracy guaranteed.

Superior quality

Analytical, with rider

Material: Stainless Steel

Application: Lab use.