

## Riders (5mg to 10mg)



**Product Categories:** [Measurements](#), [Mechanics](#), [Physics Equipments](#)

**Product Tags:** [chemistry lab instruments](#), [physics lab instruments](#)

**Product Page:** <https://www.labappara.com/product/riders-5mg-10mg/>

### Product Description

Riders (5mg to 10mg) - Physics, Chemistry Equipment

A Rider is a small piece of aluminum, from which various small weights can be attached to a balance to determine the weight of small, light objects.

#### **Features:**

Accuracy guaranteed.

Superior quality

*Analytical*, with rider

Material: Stainless Steel

**Application:** Lab use.