

## Semi Conductor Diode As A Half/ Full Wave Rectifier



**Product Categories:** [Educational equipments](#), [Electronics](#)

**Product Tags:** [electronics](#)

**Product Page:**

<https://www.labappara.com/product/semi-conductor-diode-half-full-wave-rectifier/>

### Product Description

Semi Conductor Diode As A Half/ Full Wave Rectifier - Electronics Laboratory Equipment

Unit consists of two D.C meters.

Conductor two diode & two load resistance.

Capacitor one choke & load resistance with step on transfer.