

Semi Fowler Bed



Product Categories: [Hospital Equipment & Furniture](#)

Product Tags: [hospital](#), [medical](#)

Product Page: <https://www.labappara.com/product/semi-flower-bed/>

Product Description

Semi Fowler Bed - Hospital Equipment & Furniture

Mechanically operated two sectional, uniformly perforated CRCA sheet top beds.

ABS head and foot bows / panels

Strong "H" Shaped legs on rubber stumps.

One Fix lever handles at foot end to obtain head raise upto 80deg.

Overall Size - 208L X 92W X 60H