

Senior Rotary Microtome



Product Categories: [Biology Lab](#), [Educational equipments](#), [Microtome](#)

Product Tags: [biology lab equipment](#), [science lab equipment](#)

Product Page: <https://www.labappara.com/product/senior-rotary-microtome/>

Product Description

Senior Rotary Microtome - - A microtome is a tool used to cut extremely thin slices of material, known as sections. - These Important in science, microtomes are used in microscopy - Used for the preparation of samples for observation under transmitted light or electron Microtomes use steel, glass, or diamond blades - Depends upon the specimen being sliced and the desired thickness of the sections being cut. - Steel blades are used to prepare sections of animal or plant tissues for light microscopy history. **Features:** - (Spencer 820 type) latest model, with one razor 120mm and honing stone. - Feed range 1-50 micron in steps of one micron with automatic feed release device. - Optimum performance - Less maintenance