

## Single Nose Microscope



**Product Categories:** [Microscope](#), [Microscopes](#)

**Product Tags:** [biology lab equipment](#)

**Product Page:** <https://www.labappara.com/product/single-nose-microscope-2/>

### Product Description

Single Nose Microscope - Laboratory Equipment

With one eyepiece 10 x and one objective 10 x

Total Magnification: 100 x Plano concave mirror on gamble mount for illumination, Gives Clear Images, Easy To Operate, Ideal for School & Kids. Well packed in thermometer box.

#### **Description:**

Sand casted body Inclination 90° on Horse shoe base. Focused by rack and pinion. Total Height 280 mm & Mechanical Tube Length 160 mm. Fixed square stage size 100 mm x 100 mm with stage clips. Nicely packed in Thermocol Box for packing eye piece.

## **Features:**

Smooth functionality

Long working life

Minimum maintenance

## **Optical Combination:**

Eyepiece 10 x

objective 10 x

Magnification: 100 x