

Slotted Weight (Brass)



Product Categories: [Measurements](#), [Mechanics](#), [Physics Equipments](#)

Product Tags: [physics lab instruments](#)

Product Page: <https://www.labappara.com/product/slotted-weight-brass/>

Product Description

Slotted Weight (Brass) - Physics Equipment

Set of 5 weights including hanger C.P polished.

These slotted weights are precisely turned from brass rods, carefully adjusted.

These slotted weights are supplied either in chrome plated or in brass polish finish.

Specifications:

Construction : With adjusting cavity

Material : Brass, Stainless Steel & Mild Steel

Shape : Flat Cylindrical with Slot

Uses:

Used for weighing raw materials while experiments.

Features:

Fine finish

Rust resistant coating

Durability

Capacity -

- | | | |
|----|---------|--------|
| a) | 5gmx5 | 25gms |
| b) | 10gmx5 | 50gms |
| c) | 20gmx5 | 100gms |
| d) | 50gmx5 | 250gms |
| e) | 100gmx5 | 500gms |