

Sodium Vapor Lamp Imported



Product Categories: [Educational equipments](#), [Electronics](#)

Product Tags: [optics](#)

Product Page: <https://www.labappara.com/product/sodium-vapour-lamp-imported/>

Product Description

Sodium Vapor Lamp Imported - Scientific Equipment

A sodium-vapor lamp is a gas-discharge lamp that uses sodium in an excited state to produce light at a characteristic wavelength near 589 nm.

Features:

Safe designs

Superior quality

Highly reliable

Reasonable prices

Available in watts :

(a) 35 watts

(b) 55 watts