

## Sodium Vapor Lamp Imported



**Product Categories:** [Educational equipments](#), [Electronics](#)

**Product Tags:** [optics](#)

**Product Page:** <https://www.labappara.com/product/sodium-vapour-lamp-imported/>

### Product Description

#### Sodium Vapor Lamp Imported - Scientific Equipment

A sodium-vapor lamp is a gas-discharge lamp that uses sodium in an excited state to produce light at a characteristic wavelength near 589 nm.

#### **Features:**

Safe designs

Superior quality

Highly reliable

Reasonable prices

Available in watts :

(a) 35 watts

(b) 55 watts