Soil Thermometer (brass)



Product Categories: Chemistry Equipments, Thermometer & Hygrometer

Product Tags: chemistry lab instruments, physics lab instruments

Product Page: https://www.labappara.com/product/soil-thermometer-brass/

Product Description

Soil Thermometer (brass)- Scientific Equipment

Range 20-180 F or 0-100 in body with brass cone of bottom cone 100 mm long. These thermometers are used for reading temperature of earth at depths from surface down to 1 m.

Soil thermometers have the benefit that measurements can be taken without

disturbing the thermometer.

Some of the unique features of these are given below:

Precise readings

Compact design

Quick Result

Helpful in maintaining soil temperature

Encased in corrosion free aluminum

Reads both Fahrenheit and Celsius aluminum carded

Reasonable