

Soxhlet Fat Extraction



Product Categories: [Extraction Apparatus Borosilicate Glass](#), [Laboratory Glassware](#)

Product Tags: [glassware](#), [Labappara Equipment](#)

Product Page: <https://www.labappara.com/product/soxhlet-fat-extraction/>

Product Description

Soxhlet Fat Extraction - Laboratory Glassware Equipment

Soxhlet Fat Extraction consists of R.B. Flask, Aspirator type extractor with fitting, and double surface condenser.

Fat Analyzer tests fat content according to soxhle extraction principle

It features extraction, reagent recovery, cooling system and temperature controlling.

Characteristic:

1. Small volume, fast heating speed, even heating.
2. Extraction bottle and apparatus are made of glass to avoid other spills
3. Accurate digital temperature controlling, easy to operate.

Features:

Perfectly designed

Optimum functionality

Accuracy

Precise dimensions

Easy to operate

Resistant to chemical

Lightweight

Excellent durability

Supreme finish

Necessary accessories

Precision balance 620 g/0.001 g

Available in extraction capacity- 1000 ml, 2000 ml, 3000 ml, 4000 ml, 6000 ml.